



HIRE EQUIPMENT SAFETY GUIDE

Knee Scooter

PURPOSE

A mobility device for foot / ankle injuries for patients who cannot weight bear on the affected area



FEATURES

- Patient flexes knee, and padded platform pad supports the shin of the unusable limb
- The opposite foot makes contact with the floor or ground, providing propulsion
- Folding frame for storage and transportation
- Height adjustable molded handles
- Handbrakes can be locked for added stability
- Large 8" castors for greater maneuverability and mobility outdoors
- Basket included

INSTALLATION / SETUP

Adjust handles and knee pad to the appropriate height by unscrewing the lock hubs and sliding the handles & knee pad up or down. The height should ideally be set by your therapist, but as a guide your grounded leg should be straight with your supported knee at 90 degrees.



MAINTENANCE

The knee walker should be checked monthly for signs of rust, cracks, break performance and wheel wear and tear. If any of these are present, stop use and contact your supplier.

This product should be cleaned regularly with warm soapy water unless otherwise stated by the manufacturer.

SAFETY PRECAUTIONS

- Use two hands when using the knee walker to avoid tipping.
- Ensure brakes are locked on before placing weight on the walker.
- Do not stand on this equipment.
- You should only use the equipment in the environment in which it is prescribed.
- 136kg weight limit unless otherwise stated by the manufacturer.
- Ensure the lock hubs on the height adjustable handles are secure before use.