



# HIRE EQUIPMENT SAFETY GUIDE

## Bath Board

### PURPOSE

A bath board provides a comfortable and secure platform for showering over a bathtub



### FEATURES

- Adjustable removable grab rail to assist with transfers
- Independently adjustable pillars to securely fit most baths
- Non slip surface with drain holes
- Integrated soap dish holder
- SWL 150 kg

### INSTALLATION / SETUP

To fit the bath board to your bathtub; Loosen the wing nuts located underneath the bath board. Place the board on top of the bath and adjust the width of the sliding clamps to fit the rim of the bathtub.

To use, sit on the side of the bath board that is closest to the edge of the bath with your feet on the bathroom floor. Shift or slide your bottom backwards on the board until you reach the middle. Turning slowly, hold onto the handle on the bath board and lift each leg separately over the edge of the bath, then position yourself in the centre of the bath board facing the shower. Remain seated. To exit the bath, reverse the steps above.

### MAINTENANCE

Wipe with a mild soap or detergent regularly to keep bath board clean. Holes in the seat allow water to drain quickly so sitting water will not collect.

### SAFETY PRECAUTIONS

- Supervision or assistance is recommended for people with limited mobility
- Make sure the rail is on the wall or nontransfer side of the bath
- This equipment is for sitting only
- Do not stand on this equipment
- You should only use the equipment in the environment in which it is prescribed